

Freia, (Frejya, Freyia, Frøya, Frøjya, and Freja) in Old Norse the "Lady", one of the Vanir gods, rules over the heavenly afterlife field Fólkvangr and there receives half of those that die in battle.

FREIA – a reflectometer for kinetics and liquid surfaces

H.Wacklin (ESS) & Anette Vickery (NBI)

IKON5

Lund, 2013-09-25 & 26



Freyja is also the goddess of love, beauty, fertility, gold, witchcraft, war, death...and cats



**EUROPEAN
SPALLATION
SOURCE**

Fast Reflectometer for Extended Interfacial Analysis

- Science Drivers
 - A broad spectrum of science:
 - thin film deposition/growth, structure and phase behavior
 - self-assembly and interactions of soft matter at interfaces
 - structure and interactions in biological materials and membranes
 - (thin film magnetism, interfacial effects)
 - gas-liquid, liquid-liquid, solid-liquid, solid-gas interfaces
 - (towards 2D-3D thin film structure)
- Technical Challenges and Opportunities
 - optimisation of flux and dynamic q-range for time-resolved studies
 - a range of sample sizes ($l = 5\text{-}40\text{mm}$)
 - flexible resolution options from $d\lambda/\lambda = 2\text{-}20\%$
 - horizontal sample geometry for free liquid surfaces (~25% users, up to 55% on optimised instruments e.g. FIGARO).

- Time-resolved reflectometry for interfacial kinetics primarily in
 - self-assembly of surfactants, polymers and proteins
 - rearrangement processes in thin films: e.g. lipid flip-flop, polymer interdiffusion, annealing effects in e.g. photovoltaic films
 - solvent penetration and drying/exchange effects
 - switchable materials (chemically, electrically, light driven)
 - surface reactions on the colloidal time scale; e.g. enzyme catalysis, receptor-ligand binding, drug-target interactions, surface functionalisation etc. etc. etc.
- Deuteration give unique contrast not available using other techniques (e.g. AFM, fluorescence methods, QCM-D, SPR, DPI...)
- Requirement: time-resolution in the range of ms-s-min with a dynamic q-range that allows structural analysis - unique to tof-reflectometry (cf XRR and ellipsometry, which measure a single q-point at a time).
- Challenge: typically need to reconfigure instrument 2-3 times for different angles of incidence to measure whole Q-range
- thin films on liquid surfaces can easily be disturbed by mechanical motion

The FREIA Solution

- Ven 2009: two guides?
- Lars van Moos, K.Lefmann 2011 - elliptical guide

- Deliver a broad angular divergence focused onto the sample surface
 - no need to move sample
 - angle selection by slits




- Inclined guide to access horizontal surfaces
- Angular range for Q -range for liquids ($0.005 < Q < 0.4 \text{\AA}^{-1}$) ($bkg \sim 10^{-6}$)
- Low resolution mode for high flux measurements of fast kinetics
- Higher resolution mode for thicker structures/stationary measurements
- Higher q -values can be accessed by rotation of sample


Instrument layout

- 25m instrument, elliptical guide inclined at -2° , transported divergence 4°
- $m=6$ mirror (-1°) deflects beam into guide, 4m bender ROC=56m (tbc)
- sample at 22m from moderator, 2m collimation, detector distance 1-3m

side view





top view





Bandwidth, resolution, etc.

- Aim to use natural wavelength resolution and maximum bandwidth for high flux
- e.g $2.5\text{\AA} < \lambda < 12\text{\AA}$, $d\lambda/\lambda = 3.8\text{-}18.3\%$




- 3 wavelength choppers at 6.5, 10m and 15m, 14Hz/7Hz
- 2.86ms
- Counter-rotating pairs with variable opening.




Bandwidth, resolution, etc.

- WFM challenging for 9.5Å bandwidth and 5cm guide, but possible (just):
- 7-fold WFM with $d\lambda/\lambda = 2\%$




Bandwidth, resolution, etc.

- WFM challenging for 9.5Å bandwidth and 5cm guide, but possible (just):
- 7-fold WFM with $d\lambda/\lambda = 2\%$




Simulations (without WFM)


With chopper gaps
horizontal monitor




With chopper gaps
horizontal monitor



With chopper gaps





With chopper gaps
horizontal monitor




flux on 4cm x 4cm horizontal sample (collimated)

footprint as function of wavelength and angle (uncollimated)




Simulations (without WFM)


With chopper gaps
horizontal monitor




With chopper gaps
horizontal monitor



With chopper gaps




With chopper gaps
horizontal monitor




flux on 4cm x 4cm horizontal sample (collimated)

14 Hz, $\Delta\theta/\theta=4\%$





footprint as function of wavelength and angle (uncollimated)



Collimation


- a collimating guide (2m) to clean up off-specular scattering from the guide
 - with three pairs of slits that are pre-positioned, possibility to open/close record full q-range with high time-resolution



- beam selection on pulse time scale to record full q-range 'simultaneously'.

Collimation

- a collimating guide (2m) to clean up off-specular scattering from the guide
 - with three pairs of slits that are pre-positioned, possibility to open/close with pulse-sequence and record full q-range with high time-resolution




Collimation

- Not much space available - concept for triple slit package:

3 pairs of
horizontal slit
pairs, each
movable to cover
full theta range:

a common pair of
vertical blades to
limit beam width;

slit pair position in
beam direction
should not vary
more than 40mm




Collimation

- Not much space available: use piezo-driven positioners with nm-precision


3 pairs of horizontal slit pairs, each movable to cover full theta range:

a common pair of vertical blades to limit beam width;



slit pair position in beam direction should not vary more than 40mm



SUPERACO	ATTOCORE	
SLC-17xx-S	A1P x 101xx/PES	CLOSED LOOP
X	X	TRAVEL
21, 31, 63 mm	5, 12, 30 mm	RESOL.
50 μm	50 μm	STEP SIZE
±10 μm/10 μm	2 μm	REPEATABLE
±1 μm/travel	1% full travel	ACCURACY
10 mm/s	3 mm/s	VELOCITY




piezo devices
24mm x 24mm x 11mm




Summary

- FREIA will be a versatile reflectometer for soft matter, optimised for kinetic studies on thin films
- no existing reflectometer operates without moving the sample
- Varying angle with pulse sequence has potential to allow a very broad dynamic q-range (needs development and prototyping).



enzyme kinetics slowed down by 10x in D₂O – today we can not measure this in H₂O and avoid isotope effects...



Acknowledgements



Anette Vickery
Simulations



Zvonko Lazic
Slit concept



Thomas Gahl
Motion control for slits



Ken Andersen



Iain Sutton
Choppers



Phil Bentley
Optics



Richard Hall-Wilton
Detectors

Freia, (Frejya, Freyia, Frøya, Frøjya, and Freja) in Old Norse the "Lady", one of the Vanir gods, rules over the heavenly afterlife field Fólkvangr and there receives half of those that die in battle.

Thank you!

Questions?



Freyja is also the goddess of love, beauty, fertility, gold, witchcraft, war, death...and cats



**EUROPEAN
SPALLATION
SOURCE**

Thor, the god of thunder, lightning, storms, oak trees, and strength, is one of the once rival god tribe of Æsir

