

MIRACLES STAP Report

21-22/10/2025

BENEDETTA ROSI (ESS)
On behalf of the MIRACLES team

Agenda



- 1 Guides
- 2 Choppers
- 3 Analyzer
- 4 Radial Collimator
- 5 Detector
- 6 Cave
- 7 Experimental Station
- 8 Common Projects

Overview

MIRACLES: time-of-flight backscattering spectrometer



QUICK FACTS


Sector & Beamport West (W5) Length (source to sample) 162.5 m Q-Range $0.2-2~\text{Å}^{-1}$ Wavelength Band ($\Delta\lambda$) ~ 1.5 Å (~ $\pm 500~\mu eV$) Wavelength Range 2 – 20 Å E-Transfer Range (E=2.08 meV) -2 to +20 meV

HIGH RESOLUTION QENS MODE

Energy Resolution (QENS) $\delta(\hbar\omega) \sim 2.5 \ \mu eV$ Pulse duration $\Delta t = 70 \ \mu s$

HIGH FLUX QENS MODE


Energy Resolution (QENS) $\delta(\hbar\omega) \sim 45 \mu eV$ **Pulse duration** $\Delta t = 2.4 ms$




Overview

Miracles Critical Paths





Guides




In-bunker guide system


ess

Completed in 2024


Pending: installation and integration of PWD-PS chopper module




Out-of-bunker guide system













Installations

- Thermal shutter (March 2025)
- E02 guide supports (June 2025)
- E02 guides (March 2025, June 2025, October 2025)








Swiss Neutronics

Out-of-bunker guide system




SwissNeutronics





- Integration with thermal shutter and FO chopper
- Reported instability
 - ~1 mm along x axis for guides casing
 - ~1 mm along z axis for the chopper
- Further check will be done after completion of the guide section

Choppers




Bandwidth/frame overlap (FO) chopper

Slow (Al) chopper $\alpha = 105^{\circ}$

- Installed in June 2025
- Integrated with the guide system





10


PWD-PS chopper module


Pulse width definition (PWD) chopper

- Fast (CFRP) double-disk chopper $v = 14 \text{ Hz} \dots 308 \text{ Hz}$
- $\alpha = 12^{\circ}$

Pulse shape (PS)

- Slow (Al) single-disk chopper
- $\alpha = 15^{\circ}$






Status:


• PWD-PS module manufactured and passed FAT tests








Check of performances using McStas code



Support of DMSC (Mads, Christian)

Analyser



Analyser system

Double-spherical monochromator R=2.5 m


Horizontal acceptance: 155.52°


vertical acceptance: ±22°

1032 Si(111) crystals ($\lambda = 6.27 \text{ Å}$)


Mechanical support consisting of:

- 12 carbon fiber plates (6 top, 6 bottom)
- 6 supports
- Positioning guiding system
- Alignment system










Analyser system


Assembling procedure

1. Coating with absorber (Gd_2O_3), thickness ~80 µm

2. Glueing (LOCTITE® EA 9466 $^{\text{TM}}$) onto analyzer panel (3 hours, pressure P = 0.2 MPa)










- upgrade of «coating station»
- upgrade of «assembly station» (1+2 new stations)










Analyser system

Verification and validation (preliminary)

Step 1: V&V of the final analyzer assembly


- Visual inspection
- Dimension/surface inspection (laser tracker)
- Test under vacuum






Step 2: SAT for the analyzer assembly

- Lifting/installation
- Positioning and alignment (laser tracker)
- Accuracy for sample-detector centering: ±0.15°




Radial Collimator



Radial oscillating collimator (ROC)


ess

- Aluminium structure framing carbon fiber blades coated with ~30 μm Gd₂O₃ layer
- 49 Slits (width $\varphi = 3.24^{\circ}$), angular coverage from 6.24° to 165°
- Inner radius: 30 cm
- Outer radius: 60 cm



- Manufactured by AVS
- A second layer of Gd₂O₃ will be applied to the blades at ESS Bilbao (expected to be concluded before end of 2025)

Detector









19

Angular coverage:

- 9.5° to 165° (2θ)
- ±22° vertical acceptance






96 ³He tubes (48 doublets), \emptyset =12.7 mm, P =10 bar (LND "TwinPack")



ess

- 4 preamps (2 top detector bank + 2 bottom detector bank) (CAEN)
- 1 HV supply (CAEN)
- 1 Digitizer (CAEN)




Sample area access from the roof

Cave transitable roof

Beam shielding

Cave foundation

11 1 1 1 1 1 1 1 1 1 1 1 1 1




Detectors-preamps: 48 SHV cables

- from vessel to preamps (1.5 m)
- from detectors to patch panel (~0.3 m)


Preamps-digitizer (rack): 12 RJ45 cables (~15 m)


Preamps-power supply (rack): 2 SHV power cables (~15 m)





Courtesy of O. Sidiropoulou and N. De La Rosa






Detectors-preamps: 48 SHV cables


- from vessel to preamps (1.5 m)
- from detectors to patch panel (~0.3 m)

Preamps-digitizer (rack): 12 RJ45 cables (~20 m)

Preamps-power supply (rack): 2 SHV power cables (~20 m)









Courtesy of O. Sidiropoulou and N. De La Rosa

ess

- Top preamps accessible when slab is open
- Bottom preamps accessible through the cave's maintenance door at E01 floor







Detector system Shielding

- 1 mm Cd sheet (pink)
- 5 mm B₄C-epoxy sheets (cyan)

Neutron transmittance (flux with shielding/flux without shielding)

Layout sample/ ···/detector	Point 1	Point 2	Point 3
B_4C	1.8×10^{-6}	6.0×10^{-6}	4.2×10^{-6}
B ₄ C/Cd	0.5×10^{-6}	4.2×10^{-6}	2.5×10^{-6}
$B_4C(\pm 1^\circ)$	7.8×10^{-8}	8.2×10^{-8}	9.6×10^{-8}
$B_4C/Cd/(\pm 1^\circ)$	2.2×10^{-8}	3.2×10^{-8}	2.5×10^{-8}






First steps:

- Check integration of 3he tubes with CAEN electronics
- Check of performances of each tube
 - Peak height
 - Energy resolution
 - Tests at different voltages (proportional range)
- Position sensitivity (Cd mask)

Tests performed at Lund University with Am/Be source






Courtesy of O. Sidiropoulou and N. De La Rosa



Experimental station



Cave

Cave dimensions:

- max. length: 11.690m


- max. width: 13.490m


- max. height: 8.794m

Main materials used:

- High density concrete C40/50
- Regular density concrete C40/50
- S275 JR Steel


- manufacturing phase: pre-build is expected by spring 2026
- the electrical design of the slab will be tendered





Cave






Cave

ess

Shielding (in vessel)


- 1 mm Cd sheet (pink)
- 5 mm B₄C-epoxy sheets (cyan)

To be installed end of 2025/beginning of 2026



Slab





Role: covering sample environment area during operation

- 3 high-density concrete blocks over a steel platform
- Additional layer of B₄C is integrated

The design passed Sub-TG3 review

The electrical design (including safety components e.g. switches, fence, handrails, and interfaces with PSS) will be tendered

20 mm steel structure

5 mm B₄C plate

Slab

Operation

BEAM ON:

- shutter is open
- slab is closed


When shutter is open, PSS:


- Impedes opening the fence's door
- Impedes opening the slab
- Impedes opening of the maintenance door

BEAM OFF

- shutter is closed
- slab is open

The fence's door can be opened









31

Crane



- used to transfer SE equipment from/to (1) E01 floor or (2) cave roof to the sample loading area
- technical specification document under review by ESS






Sample exposure system


Top-load

most SE equipment for MIRACLES can be hosted within the


sample chamber (inner $\emptyset = 420 \text{ mm}$)



Special cases: transition flanges



L2 inner diameter: 505 mm



L3 inner diameter: 255 mm

Sample environment

CCR/cryofurnace


Spacial considerations:

- Vertical space to host the linear exchanger
- Radial space to host the slits (and polarization cell)

Motion:

- Z linear motion (mainly to be used with the linear exchanger)
- θ rotational motion (around the Z axis)

	Range	Accuracy	Repeatability	Resolution	Speed
Z (linear)	±100 mm	±0.1 mm	±0.05 mm	0.01 mm	10 mm/s
θ (rotational)	360°	±0.1 °	±0.05°	0.01 °	2 °/s




Sample environment

Ongoing projects



Application for 2 RAC proposals

• QENS setup for in situ illumination (Forschungszentrum Jülich, Tartu Universiy, ILL, Umeå University, ESS)





 QENS setup for rapid mixing (Lund University, Tübingen University, KTH, TU Darmstadt, ESS)


Common Projects


Infrastructure CUP, CEP

- Electrical routing
- Racks layout
- Feedthroughs
- Positions of vacuum pipes, gas lines, etc.















Infrastructure

Cave feedthroughs

E01floor feedthrough



Polarization analysis capabilities




Polarizer


V-cavity, Fe/Si polarising supermirror centering at 6.27 Å


Analyzer

- Wide-angle ³He analyser cell within sample vessel
- Analyser field: Ø7m coil pair outside detector vessel, 1 mT vertical field



3He cell





Courtesy of W. T. Lee


Polarization analysis capabilities

Simulated performances



Field strength: ~ 1 mT

Field uniformity: $|\nabla \hat{b}| < 0.5 \text{ mT}$ (~ 600 hours for 3He cell)



Polarised incident beam

Flexible Braided Rigid copper Aluminium busbar Dusbar Dusbar Vessel support AISI 304


Courtesy of W. T. Lee

Thermal analysis: passive cooling (needed cooling system for HV power supply)

Beam monitors

Common Beam Monitors Project





BM1) in-bunker BM at 8.8 m, after the PWD-PS chopper module (ionization chamber)

BM2) out-of-bunker BM at 55 m, after the FO chopper (ionization chamber/bGEM)

BM3) at 162.2 m, before sample (multi-wire proportional counter)

BM4) after-sample monitor (neutron camera)

MIRACLES team













A. Conde (lead engineer)



J. Pereira



I. Aranda



I.Mazkiaran



R. Martínez







A. Simelio (in-kind)

B. Rosi





C. Beck (DMSC)





H. Bordallo



Thank you!

