




The ESTIA beamline

PRESENTED BY JOS COOPER



PSI




EUROPEAN
SPALLATION
SOURCE



Neutron Reflectometry

Instruments



~5cm

~200 μ m


Low angles mean a small slit 2 gap translates to a large beam projection onto the sample

100 μ m at 0.5 degrees = >11.5 mm footprint

Lots of lost neutrons, large footprints on samples

Neutron Reflectometry

Basics




$$n_1 \sin\theta_1 = n_2 \sin\theta_2$$

$$Q = \frac{4\pi \sin\theta}{\lambda}$$

$$n\lambda \neq 2d \sin\theta$$

Neutron Reflectometry





Neutron Reflectometry

Instruments


Time of Flight

$$Q = \frac{4\pi \sin\theta}{\lambda}$$

Monochromatic


A range of wavelengths at a fixed angle:

$$[Q_{min}, Q_{max}] = \frac{4\pi \sin\theta}{[\lambda_{max}, \lambda_{min}]}$$



A range of measured angles at a fixed wavelength:

$$[Q_{min}, Q_{max}] = \frac{4\pi \sin[\theta_{min}, \theta_{max}]}{\lambda}$$





Selene Guide Concept

Instruments


Selene guide on a ToF source:

$$Q = \frac{4\pi \sin\theta}{\lambda}$$

$$[Q_{min}, Q_{max}] = \frac{4\pi \sin[\theta_{min}, \theta_{max}]}{[\lambda_{max}, \lambda_{min}]}$$



Selene Guide Concept




Eur. Phys. J. Appl. Phys. (2012) 58: 11001
<https://doi.org/10.1051/epjap/2012110295>

Focusing specular neutron reflectometry for small samples


J. Stahn^{1a}, U. Filges² and T. Panzner²

Selene Guide Concept

Coma Aberration



The ESTIA Instrument




ESTIA

Science Case

- Magnetic systems
 - Small samples with polarisation
- *In operando* measurements
 - High flux
- High throughput measurements of e.g. biological samples
 - Low backgrounds

A reflectometer generally wants most of these features but having Selene guides and polarisation built into its heart allows ESTIA to focus on its abilities



Concept of Operations for Polarized Reflectometer (*Estia*)



Version 1.3

	Name	Affiliation
Authors	Artur Glavic, Sven Schütz	PSI
Reviewers	Jorgen Andersson Phil Bentley Thomas Gahl Francesco Piscitelli Hanna Wacklin	ESS ESS ESS ESS ESS
TG2 Review Committee	Shane Kennedy (Chair) Ken Andersen John Ankner Robert Cubitt Sylvain Désert Richard Hall-Wilton Arno Hiess Gábor László Iain Sutton	ESS ESS STAP STAP LLB ESS ESS ESS ESS
Approver	Andreas Shreyer	ESS

ESTIA



Small Samples

- The collimation on ESTIA is done at the “virtual source” – the point where we decide what neutrons to bring to the sample point (almost 30m from detector)
- 3D slit defines focus size - mirrors beam distribution on the sample point



ESTIA

Polarisation



Efficient polarization analysis for focusing neutron instruments

Jochen Stahn¹ and Artur Glavic¹

Published under licence by IOP Publishing Ltd

[Journal of Physics: Conference Series, Volume 862, International Conference on Polarised Neutrons for Condensed Matter Investigations \(PNCMI 2016\) 4–7 July 2016, Munich \(Freising\), Germany](#)

Citation Jochen Stahn and Artur Glavic 2017 *J. Phys.: Conf. Ser.* **862** 012007

DOI 10.1088/1742-6596/862/1/012007

ESTIA

Polarisation



Two transmission polarisers in series

- Improved polarisation

Part logarithmic spiral and part linear

- Linear covers rays created by coma which are almost parallel
- Pointing backwards to make device shorter
- Both sides coated with $m=5$ Fe/Si


Device is almost 1m long




ESTIA

High flux

- Divergent beam covers much larger Q space resulting in equivalent flux gain
- Moderator design leads to flux gain over existing spallation sources
- Accelerator energy leads to flux gain over existing sources



~10x flux



~10x flux





~5x flux

ESTIA

Low Backgrounds

- Elimination of line of sight three times from moderator/target
- Lots of shielding around the beam
- Enabling measurements solely of sample backgrounds













ESTIA

Science Case

- Magnetic systems
 - Small samples with polarisation ✓
- *In operando* measurements
 - High flux ✓
- High throughput measurements of e.g. biological samples
 - Low backgrounds ✓











ESTIA Science

- Spin coated samples
 - Just measure the best spots!
- Freestanding oxide membranes
 - Only measure the membranes
- Working devices
- Atmospheric chemistry
- Surfactants – soaps/lubricants
- Lipid bilayers -biology
- Composition exploration studies
 - Co-sputtering or temperature gradients
- Li ion batteries
- Complex magnetic heterostructures
- Off specular and GI-SANS



ESTIA Science

- Spin coated samples
 - Just measure the best spots!
- Freestanding oxide membranes
 - Only measure the membranes
- Working devices
- Atmospheric chemistry
- Surfactants – soaps/lubricants
- Lipid bilayers -biology
- Composition exploration studies
 - Co-sputtering or temperature gradients
- Li ion batteries
- Complex magnetic heterostructures
- Off specular and GI-SANS





**EUROPEAN
SPALLATION
SOURCE**